

Product Datasheet

Adenovirus Antibody J108

Product Name	Adenovirus Antibody
Catalog Number	J108
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified Adenovirus 5
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, purified IgB Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Does not cross-react with Parainfluenza viruses 1, 2, & 3, Influenza A & B viruses, RSV, human Metapneumovirus, and Rhinovirus
Specificity	Hexon
Preservatives	0.1% Sodium Azide