

Product Datasheet

Rubella Virus Antibody J109

Product Name	Rubella Virus Antibody
Catalog Number	J109
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Rubella virus strain HPV77
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Virions
Preservatives	0.1% Sodium Azide