

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody J111

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	J111
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	HIV-1 strain IIIB p24
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Isotype	IgG1
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Does not cross-react with human T or B cells
Specificity	p24
Preservatives	0.1% Sodium Azide