

# Product Datasheet

## Campylobacter jejuni Antibody J121

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Product Name        | Campylobacter jejuni Antibody                                                               |
| Catalog Number      | J121                                                                                        |
| Species/Host        | Rabbit                                                                                      |
| Tested Applications | ELISA, WB, IFA                                                                              |
| Immunogen           | C. jejuni ATCC strain 29428                                                                 |
| Brand               | EastCoast Bio                                                                               |
| Form/Appearance     | Purified Polyclonal Antibody, IgG Fraction                                                  |
| Storage             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.                   |
| Note                | Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available. |
| Clonality           | Polyclonal                                                                                  |
| Purity              | >95%                                                                                        |
| Buffer              | 10 mM Phosphate Buffered Saline pH 7.2                                                      |
| Specificity         | All Antigens                                                                                |
| Preservatives       | 0.1% Sodium Azide                                                                           |