

Product Datasheet

Cryptosporidium parvum Antibody J122

Product Name	Cryptosporidium parvum Antibody
Catalog Number	J122
Species/Host	Goat
Tested Applications	ELISA, IFA
Immunogen	Oocysts Purified from bovine feces
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross reacts with a large rod-shaped microbe isolated from bovine feces
Specificity	Oocysts
Preservatives	0.1% Sodium Azide