

Product Datasheet

Diphtheria Toxin Antibody J123

Product Name	Diphtheria Toxin Antibody
Catalog Number	J123
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified diphtheria toxoid
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, purified IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Purified Diphtheria Toxin Or Toxoid
Preservatives	0.1% Sodium Azide