

Product Datasheet

Epstein-Barr Virus (EBV- p18) Antibody J125

Product Name	Epstein-Barr Virus (EBV- p18) Antibody
Catalog Number	J125
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	EBV-VCA recombinant p18
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	p18
Preservatives	0.1% Sodium Azide