

Product Datasheet

Epstein-Barr Virus (EBV-EBNA-1) Antibody J126

Product Name	Epstein-Barr Virus (EBV-EBNA-1) Antibody
Catalog Number	J126
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified recombinant EBNA-1
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Isotype	IgG
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	EBNA-1
Preservatives	0.1% Sodium Azide