

Product Datasheet

Giardia lamblia (G. lamblia) Antibody J127

Product Name	Giardia lamblia (G. lamblia) Antibody
Catalog Number	J127
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified G. lamblia cysts
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Cysts
Preservatives	0.1% Sodium Azide