

Product Datasheet

Hepatitis A Virus (HAV) Antibody J128

Product Name	Hepatitis A Virus (HAV) Antibody
Catalog Number	J128
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified HAV HM175
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Specificity	Virions
Preservatives	0.1% Sodium Azide