

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody J129

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	J129
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	HBV subtypes ad & ay
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Monospecific
Specificity	ad/ay
Preservatives	0.1% Sodium Azide