

Product Datasheet

Herpes Simplex Virus 1 (HSV 1) Antibody J133

Product Name	Herpes Simplex Virus 1 (HSV 1) Antibody
Catalog Number	J133
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	HSV-1 strain F
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross-reacts with HSV-2
Specificity	Type 1
Preservatives	0.1% Sodium Azide