

Product Datasheet

Herpes Simplex Virus 2 (HSV 2) Antibody J134

Product Name	Herpes Simplex Virus 2 (HSV 2) Antibody
Catalog Number	J134
Species/Host	Sheep
Tested Applications	ELISA, WB, IFA
Immunogen	HSV-2 strain G
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross-reacts with HSV-1
Specificity	Type 2
Preservatives	0.1% Sodium Azide