

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antibody J135

Product Name Human Immunodeficiency Virus 1 (HIV 1-gp41) Antibody

Catalog Number J135

Species/Host Goat

Tested Applications ELISA, WB, IFA

Immunogen Recombinant ectodomain of gp41

Brand EastCoast Bio

Form/Appearance Purified Polyclonal Antibody, IgG Fraction

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Note Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.

Clonality Polyclonal

Purity >95%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Cross Reactivity Does not cross-react with human T or B cells

Specificity gp41

Preservatives 0.1% Sodium Azide