

Product Datasheet

Influenza A Virus Antibody J136

Product Name	Influenza A Virus Antibody
Catalog Number	J136
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Influenza A Virus, strain USSR
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Does not cross-react with Influenza B virus, RSV, Parainfluenza virus 1,2,3 or Adenovirus.
Specificity	H1N1
Preservatives	0.1% Sodium Azide