

Product Datasheet

Legionella pneumophila Antibody J138

Product Name	Legionella pneumophila Antibody
Catalog Number	J138
Species/Host	Rabbit
Tested Applications	ELISA, WB, IFA
Immunogen	Legionella pneumophila
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May cross-react with related microbes
Specificity	All Antigens
Matched Pairs	HM045
Preservatives	0.1% Sodium Azide