

Product Datasheet

Mycoplasma pneumoniae Antibody J140

Product Name	Mycoplasma pneumoniae Antibody
Catalog Number	J140
Species/Host	Rabbit
Tested Applications	ELISA, WB, IFA
Immunogen	Mycoplasma pneumoniae
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	May cross-react with related microbes
Specificity	Mycoplasma pneumoniae
Preservatives	0.1% Sodium Azide