

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody J141

Product Name	Respiratory Syncytial Virus (RSV) Antibody
Catalog Number	J141
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Human RSV
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Does not cross-react with Parainfluenza viruses 1, 2, 3, Influenza A and B viruses, or Adenovirus
Specificity	All Antigens
Preservatives	0.1% Sodium Azide