

Product Datasheet

Staphylococcus aureus Antibody J144

Product Name	Staphylococcus aureus Antibody
Catalog Number	J144
Species/Host	Rabbit
Tested Applications	ELISA, WB, IFA
Immunogen	Staphylococcus aureus
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May cross-react with related microbes
Specificity	All Antigens
Preservatives	0.1% Sodium Azide