

Product Datasheet

Leptospira Antibody J147

Product Name	Leptospira Antibody
Catalog Number	J147
Species/Host	Rabbit
Tested Applications	ELISA, WB, IFA
Immunogen	Leptospira biflexa
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2
Specificity	L. biflexa
Preservatives	0.1% Sodium Azide