

# Product Datasheet

## Human IFN gamma ELISA Set 851.560.010

|                     |  |
|---------------------|--|
| Product Name        | Human IFN gamma ELISA Set  |
| Catalog Number      | 851.560.010  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 12.5-400 pg/ml   |
| Sensitivity         | 5 pg/ml  |
| Tested Applications | ELISA  |
| Alternative Names   | IFN-g, IFN-gamma   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.  |
| Specificity         | Recognizes both natural and recombinant human IFN-g  |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |