

Product Datasheet

Streptococcus Group A Antibody J148

Product Name	Streptococcus Group A Antibody
Catalog Number	J148
Species/Host	Rabbit
Tested Applications	ELISA, WB, IFA
Immunogen	Group A streptococci
Brand	EastCoast Bio
Form/Appearance	Immunoaffinity purified Polyclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross-reactivity with other Streptococcal groups
Specificity	Group A Streptococcus species
Matched Pairs	J148 (self-pairs), J110 or HM133
Preservatives	0.1% Sodium Azide