

Product Datasheet

Adenovirus Antibody J312

Product Name	Adenovirus Antibody
Catalog Number	J312
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Immunogen	Adenovirus hexon protein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Adenovirus serotypes 1,2,3,4,5,6,7,8, 9,10,11,12, 14,16,18,20,26,31,34,35,36,37,40 & 41. No cross reactivity with Influenza A, Influenza B, RSV, Parainfluenza 1, 2 and 3, Mycoplasma pneumonia, H. pylori, and mammalian cells.
Specificity	Hexon Protein
Matched Pairs	J313 (detection)
Preservatives	0.1% Sodium Azide