

Product Datasheet

Adenovirus Antibody J313

Product Name	Adenovirus Antibody
Catalog Number	J313
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Immunogen	Adenovirus hexon protein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Reacts with at least 17 Adenovirus serotypes including types 40 & 41. No cross reactivity with Influenza A, Influenza B, RSV, Parainfluenza 1, 2 and 3, Mycoplasma pneumonia, H. pylori, and mammalian cells.
Specificity	Hexon Protein
Matched Pairs	J312 (capture)
Preservatives	0.1% Sodium Azide