

Product Datasheet

Influenza A Virus Antibody J320

Product Name	Influenza A Virus Antibody
Catalog Number	J320
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Reactive with H1 through H14 isolates
Specificity	Nucleoprotein
Matched Pairs	J320 (self-pairs) or J333 (detection)
Preservatives	0.1% Sodium Azide