

Product Datasheet

Influenza B Virus Antibody J323

Product Name	Influenza B Virus Antibody
Catalog Number	J323
Species/Host	Mouse
Tested Applications	LF, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Nucleoprotein
Matched Pairs	J323 (self-pairs), J335, J336 or J337 (capture)
Preservatives	0.1% Sodium Azide