

Product Datasheet

Human IFN gamma ELISA Set 851.560.020

Product Name	Human IFN gamma ELISA Set
Catalog Number	851.560.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	12.5-400 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
Specificity	Recognizes both natural and recombinant human IFN-g
Sample Types	Serum, Plasma, Cell culture supernatant