

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody J343

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	J343
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Immunogen	Hepatitis B surface Antigen (native)
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	No reactivity with normal human serum
Specificity	HBsAg, subtypes ad & ay
Matched Pairs	J345 (detection)
Preservatives	0.1% Sodium Azide