

# Product Datasheet

## Human TNF alpha ELISA Set 851.570

|                     |  |
|---------------------|--|
| Product Name        | Human TNF alpha ELISA Set  |
| Catalog Number      | 851.570  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 25-800 pg/ml   |
| Sensitivity         | 10 pg/ml   |
| Tested Applications | ELISA  |
| Alternative Names   | TNFa   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.  |
| Specificity         | Recognizes both natural and recombinant human TNF-a  |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |