

Product Datasheet

Human TNF alpha ELISA Set 851.570.001

Product Name	Human TNF alpha ELISA Set
Catalog Number	851.570.001
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	25-800 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.
Specificity	Recognizes both natural and recombinant human TNF-a
Sample Types	Serum, Plasma, Cell culture supernatant