

Product Datasheet

Helicobacter pylori Antibody J353

Product Name	Helicobacter pylori Antibody
Catalog Number	J353
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, IFA
Immunogen	Purified Helicobacter pylori
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Not reactive with related bacteria
Specificity	H. pylori flagellin
Matched Pairs	J353, J354 or HM1107 (detection)
Preservatives	0.1% Sodium Azide