

Product Datasheet

Herpes Simplex Virus 1 & 2 (HSV 1 & 2-gD) Antibody J357

Product Name Herpes Simplex Virus 1 & 2 (HSV 1 & 2-gD) Antibody

Catalog Number J357

Species/Host Mouse

Tested Applications LF, ELISA, IFA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity HSV-1 & HSV-2 gD

Matched Pairs J357 (self-pairs)

Preservatives 0.1% Sodium Azide