

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody J358

Product Name	Respiratory Syncytial Virus (RSV) Antibody
Catalog Number	J358
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Fusion protein of RSV strains A & B
Matched Pairs	HM1175 (detection), J358 (self-pairs), J359 or J360 (capture)
Preservatives	0.1% Sodium Azide