

# Product Datasheet

## Aspergillus Antibody JP002

Product Name	Aspergillus Antibody
Catalog Number	JP002
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC
Immunogen	Mixed species Aspergillus extract
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with related microbes (antiserum is unabsorbed)
Specificity	Soluble Proteins from Common Aspergillus Species
Preservatives	0.1% Sodium Azide