

Product Datasheet

Bacillus Antibody JP003

Product Name	Bacillus Antibody
Catalog Number	JP003
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Purified spores of B. cereus and B. subtilis
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with multiple Bacillus species (antiserum is unabsorbed)
Preservatives	0.1% Sodium Azide