

## Product Datasheet

### Candida albicans Antibody JP004

<b>Product Name</b>	Candida albicans Antibody
<b>Catalog Number</b>	JP004
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	LF, ELISA, IHC, IFA
<b>Immunogen</b>	Candida albicans
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	Cross reacts with other yeasts (antiserum is unabsorbed) No cross reactivity with human serum, urine or spinal fluid
<b>Specificity</b>	Many proteins within soluble C. albicans extract
<b>Matched Pairs</b>	HM216 (detection)
<b>Preservatives</b>	0.1% Sodium Azide