

Product Datasheet

Candida albicans Antibody JP004

Product Name	Candida albicans Antibody
Catalog Number	JP004
Species/Host	Rabbit
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	Candida albicans
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross reacts with other yeasts (antiserum is unabsorbed) No cross reactivity with human serum, urine or spinal fluid
Specificity	Many proteins within soluble C. albicans extract
Matched Pairs	HM216 (detection)
Preservatives	0.1% Sodium Azide