

Product Datasheet

Heartworm, Canine Antibody JP005

Product Name	Heartworm, Canine Antibody
Catalog Number	JP005
Species/Host	Goat
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	Dirofiliaria immitis
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with related filarial Nematodes (antiserum is unabsorbed) No cross reactivity with canine serum
Specificity	Adult worm soluble and structural antigens
Preservatives	0.1% Sodium Azide