

Product Datasheet

Chlamydia trachomatis Antibody JP007

Product Name	Chlamydia trachomatis Antibody
Catalog Number	JP007
Species/Host	Goat
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	C. trachomatis serovar L2 Purified MOMP
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross reactivity with C. psittacii and C. pneumoniae No cross reactivity with HEp-2 cells and egg yolk sac
Specificity	C. Trachomatis Major Outer Membrane Protein (MOMP)
Preservatives	0.1% Sodium Azide