

Product Datasheet

Cholera toxin Antibody JP008

Product Name	Cholera toxin Antibody
Catalog Number	JP008
Species/Host	Rabbit
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	Purified Cholera toxin
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross reacts with E. coli heat labile toxin
Specificity	Cholera toxin β -subunit
Preservatives	0.1% Sodium Azide