

## Product Datasheet

### Human TNF alpha ELISA Set 851.570.015

<b>Product Name</b>	Human TNF alpha ELISA Set
<b>Catalog Number</b>	851.570.015
<b>Description</b>	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
<b>Range</b>	25-800 pg/ml
<b>Sensitivity</b>	10 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	TNFa
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.
<b>Specificity</b>	Recognizes both natural and recombinant human TNF-a
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant