

Product Datasheet

Enterococcus species Antibody JP010

Product Name	Enterococcus species Antibody
Catalog Number	JP010
Species/Host	Rabbit
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	Enterococcus faecium
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with related microbes (antiserum is unabsorbed)
Specificity	Reactive With Whole Cells
Preservatives	0.1% Sodium Azide