

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody JP012

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	JP012
Species/Host	Rabbit
Tested Applications	LF, ELISA, IHC, IFA
Immunogen	subtypes Ad and Ay
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross-reactivity with human serum proteins
Specificity	HBsAg
Preservatives	0.1% Sodium Azide