

Product Datasheet

Hepatitis C Virus (HCV) Antibody JP016

Product Name	Hepatitis C Virus (HCV) Antibody
Catalog Number	JP016
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Recombinant NS5a + NS5b
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	NS5a + NS5b
Preservatives	0.1% Sodium Azide