

Product Datasheet

Herpex Simplex Virus (HSV) Antibody JP017

Product Name	Herpex Simplex Virus (HSV) Antibody
Catalog Number	JP017
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	HSV strain F (human)
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	HSV-1 and HSV-2
Specificity	ICPs and late structural proteins
Preservatives	0.1% Sodium Azide