

## Product Datasheet

### Influenza A Virus Antibody JP020

<b>Product Name</b>	Influenza A Virus Antibody
<b>Catalog Number</b>	JP020
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	LF, ELISA, WB, IHC, IFA
<b>Immunogen</b>	Influenza A, strain Texas 1/77
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	No reactivity with Influenza B Virus, RSV, Parainfluenza Viruses 1-3 or Adenovirus; possible reactivity with chicken cellular proteins
<b>Specificity</b>	Purified Virions (H3N2)
<b>Preservatives</b>	0.1% Sodium Azide