

Product Datasheet

Klebsiella pneumoniae Antibody JP021

Product Name	Klebsiella pneumoniae Antibody
Catalog Number	JP021
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Klebsiella pneumoniae
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with related microbes (antiserum is unabsorbed)
Specificity	Polyvalent for Multiple Klebsiella Antigens
Preservatives	0.1% Sodium Azide