

Product Datasheet

Lipid A Antibody JP022

Product Name	Lipid A Antibody
Catalog Number	JP022
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Lipid A from E. coli
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Reactivity with many Enterobacteriaceae
Specificity	LPS
Preservatives	0.1% Sodium Azide