

Product Datasheet

Human CD62E ELISA Set 851.580.001

Product Name	Human CD62E ELISA Set
Catalog Number	851.580.001
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	1-32 ng/ml
Sensitivity	0.5 ng/ml
Tested Applications	ELISA
Alternative Names	E-selectin / ELAM-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble ELAM-1
Sample Types	Serum, Plasma, Cell culture supernatant