

Product Datasheet

Mycobacterium tuberculosis (TB) Antibody JP023

Product Name	Mycobacterium tuberculosis (TB) Antibody
Catalog Number	JP023
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Purified protein derivative (PPD) of M. tuberculosis
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Isotype	IgG
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	Reactive with several related microbes including M. avium, M. phlei, and M. parafortuitum (antiserum is unabsorbed)
Specificity	M. tuberculosis
Preservatives	0.1% Sodium Azide